

FULLIFE HEALTHCARE PVT.LTD.




QUALITY CONTROL DEPARTMENT

FINISHED PRODUCT CERTIFICATE OF ANALYSIS



Product Name	Fast & Up Man Extend Delay Effervescent Tablet Fizzy Peach	A.R. No.	FKH/2526/QCP/NFG/01182
FG Code	NFG000256	Mfg. Date	12/2025
Batch No.	NM0925003	Exp. Date	05/2027

Sr. No	Test	Specification	Results
1.	Description	Slight peach color, mottled, round, flat, effervescent tablet plain on both sides.	Slight peach color, mottled, round, flat, effervescent tablet plain on both sides.
2.	Odour	Should give characteristics peachy apple odour.	It give characteristics peachy apple odour.
3.	Solution Taste	Should have sweet and peachy apple taste.	Solution have sweet and peachy apple taste.
4.	Appearance of solution	Should be peach color fizzy solution	peach color fizzy solution
5.	Average Weight	4200 mg \pm 3 % (Between 4074 mg to 4326 mg)	4156.6 mg
6.	Uniformity of Weight	Maximum 2 tablet / 20 tablets =Not less than 3990 mg or not more than 4410 mg Not less than 3780 mg or not more than 4620 mg	Minimum=4092.8 mg and Maximum=4218.8 mg None deviates from this limit Not Less than 3780 mg and Not more than 4620 mg
7.	Hardness	Not Less Than 60 N	149 N
8.	Loss On Drying	Not more than 1.0 %	0.40 %
9.	pH of Solution	Between 3.0 to 6.0	5.74
10.	Fizz Time	Not More Than 5 Minutes	03 Minutes 56 seconds

Signatories	Prepared By	Reviewed By	Approved By
Sign & Date	 05/01/2026	 05/01/2026	 05/01/2026

FULLIFE HEALTHCARE PVT.LTD.

QUALITY CONTROL DEPARTMENT

FINISHED PRODUCT CERTIFICATE OF ANALYSIS



Product Name	Fast & Up Man Extend Delay Effervescent Tablet Fizzy Peach	A.R. No.	FKH/2526/QCP/NFG/01182
FG Code	NFG000256	Mfg. Date	12/2025
Batch No.	NM0925003	Exp. Date	05/2027
Sr. No	Test	Specification	Test Method
11.	Identification	<p>L- Arginine: Retention time of test solution matched with retention time of reference solution.</p> <p>Black Pepper : Retention time of test solution matched with retention time of reference solution.</p>	<p>L- Arginine: Retention time of test solution matched with retention time of reference solution.</p> <p>Black Pepper : Retention time of test solution matched with retention time of reference solution.</p>
12.	Assay		
	L-Arginine -175 mg	Between 90.0 % to 110.0 %	102.8 %
	Black pepper -5 mg	Between 90.0 % to 110.0 %	102.0 %
	Magnesium-50 mg	Between 90.0 % to 110.0 %	95.2 %
	Zinc-8 mg	Between 90.0 % to 110.0 %	97.13 %
13.	Heavy Metals		
	Lead	NMT 2.5 ppm	0.047 ppm
	Mercury	NMT 0.25 ppm	BLQ
	Arsenic	NMT 1.1 ppm	0.023 ppm
	Cadmium	NMT 1.5 ppm	BLQ
14.	Microbial Contaminants		
	Total Aerobic Microbial Count	NMT 1000 cfu/g	20 cfu/g
	Total Yeast / Mould Count	NMT 100 cfu/g	10 cfu/g
	E. Coli	Should be absent/g	Absent
	S. aureus	Should be absent/g	Absent
	Pseudomonas aeruginosa	Should be absent/g	Absent
	Salmonella	Should be absent/10g	Absent

Remark: Test, Heavy metals test value incorporated from NABL Reliable Lab report references No.: RAL/FV251216033“BLQ”- Below Limit of Quantification.

Signatories	Prepared By	Reviewed By	Approved By
Sign & Date	 05/01/2026	 05/01/2026	 05/01/2026